



Lorne Buffered Saline Tablets

The Department of Health and Social Security Notice HN (Hazard) (83) 6⁽¹⁾ and other published studies^(2,3) warn against the danger of using certain saline preparations in the antiglobulin test, particularly those intended for irrigation and intravenous use. Solutions of unexpectedly low pH and those autoclaved and stored in plastic containers could severely compromise the sensitivity of the antiglobulin test, resulting in failure to detect clinically significant antibodies. This could prove fatal in certain cases.

Re-sealable pack

Lorne Buffered Saline Tablets

Each Tablet in Solution Contains 0.9% NaCl in 1 Litre of Deionised Water at pH 6.90 ± 0.10

Contains no preservative

Store tablets in a cool dry place Directions for use

Dissolve each tablet in 1 Litre of deionised Water Dissolution time 20 minutes (approximately) Store at 2-8°C when in solution

LOT NO: EXPIRY:

Re-sealable pack

Lorne Laboratories have introduced tablets containing a phosphate buffer, which produces an isotonic saline at pH 6.9±0.1 suitable both as a washing solution and also as a medium for cell suspension. The daily use of Lorne Phosphate Buffered Saline Tablets provides a safe, standardised saline solution for transfusion serology. The tablets have been quality controlled for pH, conductivity and osmolality. Once reconstituted in deionised water, the PBS solution is stable for 6 months at 2–8°C provided that the solution is stored in a clean container which is properly capped. If the PBS solution is used regularly, it may be stored for up to 7 weeks at room temperature.

 The Department of Health & Social Security, Health Services Management, Antiglobulin test false negative results. HN (Hazard) (83) 6. Nov 25. 1983.

AVAILABLE PACK SIZES

Pack Size	Catalogue Number
25 tablets	490025

Pack Size	Catalogue Number
250 tablets	490250

Lorne Labs Limited, Unit 1 Cutbush Park Industrial Estate, Danehill, Lower Earley, Berkshire, RG6 4UT, UK

Document reference number: Pl490

Document issue number: 8/02/2018

Page 1 of 1