


Flatbed Agitators

The flatbed agitator keeps your thrombocytes in motion



Agitator 96



Agitator 60



Agitator 48



Agitator 24

Storage of thrombocytes – standstill or motion

- Upper section moves over 30mm at ~60 strokes at 50Hz / ~72 strokes at 60Hz per minute
- Every shelf of Agi 24, 48, 60 has a capacity of 6 bags
- Every shelf of Agi 96 seizes 12 bags
- **Built-in stopper** prevents accidental drawing out
- Designed in **stainless steel**
- Each drawer can be **operated independently** from other drawers
- Designed for **silent and reliable long time use**

Optional

No-motion alarm

- Acoustic warn signal
- Via D-sub 9 connector relais contact

Technical data	Flatbed Agitator 24	Flatbed Agitator 48	Flatbed Agitator 60	Flatbed Agitator 96
Capacity (bags):	24	48	60	96
Power supply:	220 V / 50 Hz / 50	220 V / 50 Hz / 50	220 V / 50 Hz / 50	220 V / 50 Hz / 50
Weight (kg):	20	26,5	30	38,5
Dimensions (cm):	37,7 x 44,8 x 28	37,7 x 44,8 x 40	37,7 x 44,8 x 46,3	88,9 x 35,6 x 40,5

Lmb Technologie GmbH

Mösistr. 17
D-85445 Schwaig

Tel: +49-8122-88096-0
Fax: +49-8122-88096-60

www.lmb.de
info@lmb.de