

Climatic chamber CLIMAX

Now also available as
medical ———
————— **device**



CLIMAX 100



CLIMAX 120

Now with complete TEMPERATURE TRACKING

Climax climatic chambers keep together with fitting Agitators your thrombocytes in motion and at constant 22°C. The compressor of Climax series works extremely efficient to keep energy consumption at a very low level. Strong air recirculation and connected quick temperature adaption also account for the low energy consumption.

The included temperature alarm system works independently from the temperature controller (double temperature control).

Options

- Web light server allows access to the device from a computer and multiconnection of refrigerators
- Remote panel as option
- PT 100 extension module in combination with PT 100 sensors enables calibration of the device
- GSM Communication module (for sms communication)

Characteristics

- **Housing** made of double-side epoxy laminated steel
- **Interior:** made of stainless steel (18/19 AISI 304), exterior laminated in white
- Small **rolls** with fixing break
- Produced of **chrome-plated steel**
- **Temperature control** 22 °C +/- 1 °C.
- Electronically controlled **heating and cooling system**
- **Air circulating fan** for the fast assimilation of the target temperature
- **Isolation:** the cavity between the two metal walls is packed 5 cm with the foam material polyurethane
- **Hung-in Glass Door:** 3x isolation glass embedded in the aluminium frame
- **Framework:** lined with refrigerator seal; external door handle
- **Stained steel plating:** as footprint for Platelet agitator
- **Prepared inner wall:** for placing of an external measurement probe
- **Integrated Temperature Recording:** including separate temperature probe