

## QuickSlide Plus II Stainer



### Dimensions/Weight:

Width 30.48 cm (12.0 in.)  
Height 27.9 cm (11.0 in.)  
Depth 27.94 cm (11.0 in.)  
Weight 4.19 Kg (9.25 lb)

### Power Requirements:

Input 100-240 50/60HZ  
International Power Converter Available

### Temperature, Ambient Operation:

16-32 Degrees C (60-90 Degrees F)

### Humidity:

0 to 95% without condensation

### Operator Adjustments:

Stain/Buffer Times  
Adjustable Fill Levels  
On & Off Audible Sound Control

### Operator Entry:

Moisture Proof, Soft Touch Keypad  
Stain Resistant Injected Moulded Cabinet  
Special formulated reagent kit - 125 tests per kit in one self-contained package.

Lorne Laboratories Limited  
Unit 1 Cutbush Park Industrial Estate  
Danehill, Lower Earley  
Berkshire RG6 4UT United Kingdom

Tel: +44 (0) 118 921 2264  
Fax: +44 (0) 118 986 4518  
Email: [info@lornelabs.com](mailto:info@lornelabs.com)  
[www.lornelabs.com](http://www.lornelabs.com)



The QuickSlide Plus II stainer represents a significant advancement in automated slide staining and is a valuable addition to any laboratory, whether it is the primary stainer used by small hospitals, physicians labs, stat labs, or used as a stat/backup instrument in larger labs. The QuickSlide Plus II can also be utilized in veterinarian labs for a variety of biological specimens.

**New Feature!** In addition to the soft agitation before the Rinse is drained, QuickSlide Plus II has an added soft agitation feature that occurs before the Buffer is drained. Through extensive testing we have determined this added feature puts off a cleaner slide with the same great staining characteristics of our modified Wright Giemsa stain formula.

The unit needs no daily maintenance and only requires a tubing change once or twice annually.

### **Small Compact Size**

The QuickSlide Plus II is one of the smallest automated slide stainers available. It requires only one square foot of counter space. In addition to its small size, the QuickSlide Plus II reagent system is included in the one square foot.

### **Fast Staining**

The QuickSlide Plus II can stain a blood smear in as little as 60 seconds. In addition to the rapid staining, the QuickSlide Plus II eliminates cross contamination seen in the dip method and yields consistent results no matter who stains the slide.

### **Performs Multiple Staining**

In addition to staining peripheral blood smears, the QuickSlide Plus II has a special program that allows it to stain bone marrow samples. The QuickSlide Plus II also allows operator to adjust both stain and buffer times.

### **Features of the QuickSlide Plus II:**

- State-of-the-art microprocessor, increased speed and ease of operator use
- Operator prompt LCD display
- Soft touch color-coded keys
- Adjustable stain and buffer timing - so important for consistent lab results with multiple users
- Adjustable fill level
- “Cycles Remaining” display - computer monitors reagent and prompts when to reorder
- Injection moulded case with built-in reagent platform
- Stainless steel reagent cannulas - oxidation resistant, increase system durability
- Color-coded reagent tubing for ease of operation
- Special formulated reagent kit - 125 tests per kit in one self-contained package.