

## Quality Mixer



Quality Mixer is a new fully automatic mixer/tube stripper manufactured to simplify and standardise the requirements of QC according to Directive 2004/33/EC, Article 6 Appendix V2

Quality Mixer prevents ergonomic injuries in hand, wrist and shoulders which happen during manual processing. It strips the tube and mixes the bag automatically. It is adjustable from 3-10 stripping cycles and 3-10 mixing cycles. For example, the tube will be stripped five times and the bag mixed ten times between each strip. Quality Mixer is suitable for all blood bags on the market

Quality Mixer:

- Fully automatic mixer/tube stripper
- Simplifies and standardises QC requirements for blood collection
- Assists compliance with Directive 2004/33/EC, Article 6 Appendix V2
- Prevents ergonomic injuries in hand, wrist and shoulders during manual processing
- Strips the tube and mixes the bag automatically
- Adjustable from 3-10 stripping cycles
- Adjustable from 3-10 mixing cycles
- Suitable for all blood bags on the market

<b>Quality Mixer</b>	<b>150101-01</b>
Mains unit:	Input 24 V
Mixing cycle:	11 cycles per minute
Dimensions:	290 (L) x 253 (W) x 150 (H) mm
Weight:	5 kg
Max weight:	500 g
Power supply unit:	Input 230 VAC +/- 10%, output 24 V
CE marking:	CE marked